



PRODUCT SPECIFICATION



Product: Elastic bandage
REF: 1882



1 General description

An elastic bandage with firm edges equipped with a fastening clip.

2 Intended use

The Elastic Bandage is intended to be used for strong fixation bandaging, for support and relief bandaging, for fixation of a wound dressing, and for aftercare of distortion, dislocations, contusions, and fractures. Not intended for use directly on open wounds.

3 Technical description

Product dimensions

Bandage 4m x 8 cm

Material

Bandage 65% cotton and 35% polyamide
Package Polypropylene film.

4 Composition

Consumer package

<u>Ref#</u>	<u>Product/Material</u>
51618050	Elastic bandage

<u>Pcs</u>
1
1

Transport package

Consumer packages (pcs):	20
Dimensions (l x w x h) cm:	26 x 20,5 x 17,5
Weight (kg):	0,63

5 Labeling

On product/ consumer package

Cederroth address, telephone- and faxnumber, web-address, E-mail-address
Product name, Reference code
Lot-number, Expiry date, EAN-code, CE

On transport package

Cederroth address, telephone- and faxnumber, web-address
Product name, Reference code
Lot-number, Expiry date, EAN-code, Quantity

6 Storage conditions

The product should be stored in the original package for maximum protection until time of use.
Shelf life: 5 years from production date

Authorization: