

PRODUCT SPECIFICATION

REF 490920 CEDERROTH FIRST AID STATION



1 GENERAL PRODUCT DESCRIPTION

Description

A First Aid station containing a variety of first aid products, intended to be used in first aid situations. All products are easily available and protected against dust and moisture by transparent lids. Space for additional products is provided behind a lid with first aid instructions. Adhesive bandages (plasters) and wound cleansers are placed in dispensers for easy access and theft-proof storage. A special key is provided for refill replacement.

Medical devices in the station are CE-labelled.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product</u>	<u>Qty</u>	<u>REF</u> <u>(refill)</u>	<u>Product</u>	<u>Qty</u>	<u>REF</u> <u>(refill)</u>
First Aid station (empty)	1	-	Salvequick Wound Cleanser 20-p	1	323700
Cederroth 4-in-1	2	1910	Protection kit	1	2596
Cederroth 4-in-1 Mini	3	1911	Cederroth Burn Gel Dressing 2-p	1	901900
Salvequick Plastic Plasters (45) in dispenser	1	6036	Cederroth Eye & Wound Cleansing Spray	1	726000
Salvequick Textile Plasters (40) in dispenser	1	6444	First Aid Instruction (label)	1	-
			Key to the dispensers	1	51616620

Material

First aid station:	Green parts:	Polystyrene, ABS-plastic
	Transparent lid:	Polycarbonate
	Sealing profiles:	EPDM rubber (synthetic rubber)
First aid products:	Refer to Product specifications for each product	

Shelf life: 5 years from production date for each product in the station.

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310804909205	7310804909205	7312024909201
Content:	1	1 CU	11 SKU/layer / 4 layers
Size (w x l x h mm):	294 x 572 x 123	294 x 572 x 123	Pallet 800 x 1200 x 1326
Gross weight:	3153 g	3153 g	164 kg
Packing material:	Carton 313 g	Carton 313 g	-

