

PRODUCT SPECIFICATION

REF 6496 SALVEQUICK TEXTILE FINGER BANDAGE



1 GENERAL PRODUCT DESCRIPTION

Description

The Textile Finger bandages are packed in a cardboard refill intended for the Cederroth (Salvequick) dispenser for adhesive bandages.

The finger bandages are extra- long and thus well suited to be wrapped around fingers. The finger bandages can easily be pulled out with just one hand and applied on the wound.

The finger bandage consists of a backing, evenly spread with an adhesive coating, with an absorbent pad fixed centrally to the dressing. The adhesive is covered with a release-coated foil that consists of two parts partly overlapping each other. When the adhesive bandage is pulled out from the package, the lower cover will remain behind. Thus one half of the finger bandage can be put on the skin with one hand and after removing the rest of the cover the plaster is ready for use.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity (pcs)</u>	<u>Dimensions (mm)</u>	
Refill packages	6	30 x 88 x100	
<u>Per refill package:</u>		<u>Dressing</u>	<u>Wound Pad</u>
Finger bandage	30	120 x 20	24 x 20

Material

Film/backing:	Elastic, beige viscose/polyamide fabric with hydrophobic treatment
Adhesive	Acrylate
Wound pad	Viscose, polypropylene, polyethylene with polyethylene wound contact layer
Facing	Siliconized paper and plastic
Protective paper	Coldseal paper coated with latex-based adhesive
Wrapping	Carton

Shelf life: 5 years

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310610064969	7310806496000	-
Content:	6 refills	10 CU	11 SKU/layer x 5 layers
Size (w x l x h mm):	150 x 100 x 90	155 x 513 x 190	Pallet 800 x 1200 x 1100
Gross weight:	200 g	2237 g	148 kg
Packing material:	Carton 28 g	Corrugated board 237 g	-



Version: 05
Issued: 2017-01-13